

Exporting XSTD Files from i-CAT Classics

XSTD files are the native format used by Xoran software on the i-CAT Classic. You can export an XSTD file from Xoran and send it to IDT for processing. This is faster and easier than sending us the DICOM files.

Step 1

Double click the Xoran icon shown below on the desktop:




Step 2

Enter the Username and Password: - normally **ct**




Step 3

On the patient database window click the **Find** button as shown below




Step 4

Click the  next to any patient name and select a **Recon** data set as shown below. **DO NOT SELECT RAW DATA.** You do not need to load the study.



Step 5

The following dialogue box will appear:



- Browse to a convenient folder such as the **Desktop**
- Click the **Export** button

Step 6

The .xstd file will be saved to your Desktop:



Step 7 (Optional)

Use Method A or Method B or Method C or any other method to send the XSTD File to IDT.

Method A: if you have an sFTP Account

You can drag the .xstd file into the "To_IDT" folder.

Step 8A

- Find the file you made in **Step 6**
- On your Desktop, open **WinSCP**
- Drag the .xstd file you created into the **To_IDT** folder.
- When the .xstd file has finished copying over, delete it from your Desktop.

Method B: if you have registered on www.idtscans.com

You can upload the .xstd file while booking in the work that you require.

Step 8B

- Go to www.idtscans.com
- Log into your account and start booking the case in
- Select **".xstd, .inv, .spr, .sim or .s1s file"** and click on **"Upload Now"**.



The screenshot shows the IDT website interface. At the top left is the IDT logo with the tagline "making the most of your image". The main header includes the website URL "www.idtscans.com" and navigation links: "Home | About Us | Services | Support | Knowledge Base | Contact Us |". Below this is the user's name "Dr. Reynolds" and three links: "My Profile", "My Home", and "Logout". A navigation bar contains "Back" and "Next - Patient Details" buttons. The main content area is titled "How will you send the scan ?" and lists four options with radio buttons: ".zip file containing DICOM images", ".xstd, .inv, .spr, .sim or .s1s file", "DICOM images from a CD", and "IDT already has the data". Each option has a "Click Here" button next to it. The second "Click Here" button is highlighted with a red border. To the right of each "Click Here" button is a yellow information icon (a lowercase 'i'). At the bottom left, it shows "Price so far: GBP £ 0.00" and a "Refresh" button. The footer contains copyright information "©2008-17 IDT Scans", links for "Terms & Conditions" and "Privacy Policy", and the version number "V3.31 Rev 2017-12-09".

- Attach and send your file (see **Step 10** below)
- Carry on booking in the case while the file is being sent to IDT.


Method C: if you haven't registered yet

You can upload the .xstd file first, then book in the work you require later.

Step 8C

Go to **www.idtscans.com**

If you haven't registered, click 'Upload Your Scan Now' on the homepage (see below):



Step 9C

On the next page, click on "**Upload a Single File**":



Upload Files

If you are requesting work for IDT Scans to do, please click on one of the following links and follow the steps in the Online Booking wizard:

- [Reformat an Existing Scan](#)
- [Request a Radiology Report](#)
- [Implant Simulation, Treatment Planning, Guides](#)

If you want to send us a single file, please click here:

- [Upload a Single File](#)

If you want to send us multiple files, please place them in a folder. Right-click on the folder and select "Send to" and then select "Compressed (zipped) folder". Then click on the following link and select the zipped folder:

- [Upload a Zipped Folder](#)

Please visit our extensive [Knowledge Base](#) for further information, or [Contact Us](#).

Step 10

Follow the 3 steps below to attach and send your file. If you cannot see the **Browse** button please install Adobe Flash Player from <http://get.adobe.com/flashplayer>



Mail a Big File to idtscans.com

In **Secure Mode**

1 **Attach Files** *(maximum file size 2GB)*

2 **Your Email Address**

Message

it is easy to send data to
IDT Scans :)

3

Created 2015-02-01 RAR
Revised 2017-12-06 DF
Revised 2020-03-07 DF/RAR