



Downloading and Saving Simplant files that you received from IDT

The links sent to you by IDT are valid for 60 days only.

Simplant software looks for files in the C:\DentalData folder so you should download and save your Simplant files in this folder. The way to do this depends on which Web Browser you are using.

Step 1

Left click on the link (with the **Left** mouse button) to access your file:

The screenshot shows the IDT Scans website interface. At the top left is the IDT logo with the tagline "making the most of your image". In the center is the URL "www.idtscans.com". On the top right is a small image of a person's face. Below the header, there is a message: "Please use the link below to retrieve your file from IDT Scans". This is followed by instructions: "Please save this file to C:\DentalData, then open it with SIMPLANT". The filename is listed as "2017-04-07 Mn Mx JComp.sim" with a redacted area before the date. The file size is "5.3 Mb". A note states "The file will be available for 60 days". A red-bordered box highlights the download link: "Please use this link to retrieve the file: https://www.mailbigfile.com/5d19689dac7f2dd2a8d5c5af73a89b0f/listFiles.php?repro_id=3298". At the bottom, it says "To book another case please visit IDT Scans".

Step 2


Enter the password assigned to you by IDT and left click **Submit**.

The screenshot shows the IDT Scans website interface. At the top left is the IDT logo with the tagline "making the most of your image". In the center is the URL "www.idtscans.com". On the top right is a small image of a person's face. Below the header, there is a message: "Here is the file sent to you by info@idtscans.com". Below this, there is a prompt: "The sender has requested you enter a password to retrieve your file." Below the prompt is a text input field containing "*****". At the bottom is a pink "Submit" button.

If you are using Microsoft Edge:

Step 3A

Right Click (with the **Right** mouse button) on **Download** and select **Save target as ...**




Please skip to Step 4.

If you are using Firefox:

Step 3B

Right Click (with the **Right** mouse button) on **Download** and select **Save link as ...**



Please skip to Step 4.

If you are using Google Chrome:

Step 3C


Right Click (with the **Right** mouse button) on **Download** and select **Save link as ...**



All Browsers:

Step 4

Navigate to **Computer** then **Local Disk (C:)** then **DentalData** then click **Open**




Please Note: Some websites will change the filename extension “.sim” to “.zip”.

If this happens you must rename the file as follows:

Simplant Planner:	change .zip to .sim
Simplant Pro:	change .zip to .spr
Simplant OneShot:	change .zip to .s1s
Simplant Master:	change .zip to .sma


Step 5

Click **Save** to save your file in the C:\DentalData folder.



Step 6

Launch Simplant by double-clicking the Simplant icon on your Desktop.
From the **File** menu select **Open Files ...**



A dialog box will open and you will see your patient listed. Double-click on your patient to view the images.